


ABPM Tutorial – Preparing the ABPM

(Please note – ABPM Software must be installed on your computer before continuing with this guide.)

Step 1. Use the provided USB cable to connect the ABPM to your computer.

Step 2. After a short moment, you will see a small window pop up requesting your credentials. **Please enter your username and password.** Then click “New Patient”



Cardio Study - Device Connected


Email

Password

Test Completed **New Patient** Cancel

Step 3. On the next window, enter the following information:

- Patients First Name
- Patient’s Last name
- Health Card (no version code necessary).
- Awake Time (When patient usually wakes up)
- Asleep Time (When patient is usually going to bed)
- Awake Intervals (How often the monitor should take a reading while patient is awake)
- Asleep Intervals (How often the monitor should take a reading while patient is asleep)



Cardio Study - New Patient

First Name


Last Name

Health Card

	Time (24hr)	Interval (min)
Awake	<input type="text" value="8:00"/>	<input type="text" value="30"/>
Asleep	<input type="text" value="21:00"/>	<input type="text" value="60"/>

Next Cancel

Step 4. You will see a window indicating that your patient has been registered. Click "OK".



Step 5. You may now unplug the machine from your computer.